

# DENTAID

## The Oral Health Experts



VADEMECCUM

# interprox.


Cleans all spaces for complete oral hygiene

## INDICATIONS

Interdental spaces




Implants



Diastema



Orthodontic appliances



Fixed prosthesis



Furcal areas



## Quality and Hygiene

- Plastic coated wire that protects teeth, gums and dental implants
- High quality black and white Tynex® filaments for detecting plaque and bleeding
- The cap protects the brush and guarantees optimal hygiene after use

## interprox.



## interprox. plus



### Flexibility

- The head and the handle are flexible and do not have any "memory" effect allowing easier access to every interdental space
- Plastic lined neck for added protection, flexibility and anchoring of the wire

### Comfort and precision

- Anti-slip groove for precision and comfort during brushing

### Specific design

- Angled head for better access to posterior spaces
- Longer and rigid handle

# interprox.

Cleans all spaces for complete oral hygiene

## SIZES RANGE

### interprox.



### interprox. plus

## SAMPLING BOX



### interprox.gel

gel toothpaste



		Gel
Cetylpyridinium chloride	Reduces biofilm accumulation.	0.3%
Sodium fluoride	Reinforces enamel.	1450ppm
Zinc lactate	Increases the substantivity of CPC, lengthening its activity.	0.14%
Allantoin	Promotes healthy gums.	0.10%
Sizes		20ml

#### REFERENCES:

1. Schmidt B, Jentsch H. Comparison Between Two Interproximal Cleaning Procedures in Periodontitis Patients: A Six-month, Single-blind, Randomised Controlled Clinical Trial. Oral Health Prev Dent. 2015;13(3):205-11

\*ISO 16409: International standard for manual interdental brushes



Follow us on:

